


## Needleless Sets



**P/N: A2-80070**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller clamp (4), Male Luer Lock Adaptor (5), Luer Lock Tip protector (6)



**P/N: A2-80071-P**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Needleless Injection Site (5) at proximal side, Roller clamp (6), Male Luer Lock Adaptor (7), Luer Lock Tip protector (8)



**P/N: A2-80071-D**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), Slide clamp (5), Needleless Injection Site (6) at distal side, Male Luer Lock Adaptor (7), Luer Lock Tip protector (8)



**P/N: A2-80072**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (8), Roller Clamp (6), Slide clamp (7), Male Luer Lock Adaptor (9), Luer Lock Tip protector (10)

## Needleless Sets



**P/N: A2-84070-N056**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller clamp (4), Male Luer Lock Adaptor (5), Vented Cap with Filter (6)




**P/N: A2-84071-P046**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Slide clamp (5), Needleless Injection Site (6) at proximal side, Roller clamp (7), Male Luer Lock Adaptor (8), Vented Cap with Filter (9)



**P/N: A2-84071-D058**


105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), Slide clamp (5), Needleless Injection Site at distal side (6), Male Luer Lock Adaptor (7), Vented Cap with Filter (8)



**P/N: A2-84072-N048**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Slide clamp (5)&(8) Two Needleless Injection Sites, one at the proximal side (6), one at the distal side (9), Roller Clamp (7), Male Luer Lock Adaptor (10), Vented Cap with Filter (11)

# Needleless Sets




**P/N: A2-84072-N047**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (8), Roller Clamp (6), Slide clamp (7)&(11), Male Luer Lock Adaptor (9), Vented Cap with Filter (10)


## Needleless Sets



**P/N: A2-80071-P046**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Slide clamp (5), Needleless Injection Site at proximal side (6), Roller clamp (7), Male Luer Lock Adaptor (8), Luer Lock Tip protector (9)



**P/N: A2-80072-N048**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Slide clamp (5)&(8), Two Needleless Injection Sites, one at the proximal side (6), one at the distal side (9), Roller Clamp (7), Male Luer Lock Adaptor (10), Luer Lock Tip protector (11)


## Needleless Sets with Anti-Free-Flow Valve



**P/N: A2-83070**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), Male Luer Lock Adaptor (5), Male Luer Lock Anti-Free Flow Valve (6), Luer Lock Tip protector (7)



**P/N: A2-83071-P**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Needleless Injection Site (5) at proximal side, Roller Clamp (6), Male Luer Lock Adaptor (7), Male Luer Lock Anti-Free Flow Valve (8), Luer Lock Tip protector (9),



**P/N: A2-83071-D**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), Slide clamp (5), Needleless Injection Site (6) at distal side, Male Luer Lock Adaptor (7), Male Luer Lock Anti-Free Flow Valve (8), Luer Lock Tip protector (9)




**P/N: A2-83072**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (8), Roller Clamp (6), Slide clamp (7), Male Luer Lock Adaptor (9), Male Luer Lock Anti-Free Flow Valve (10), Luer Lock Tip protector (11)


## Needleless Sets with Anti-Free-Flow Valve



**P/N: A2-83071-P046**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Slide clamp (5), Needleless Injection Site at proximal side (6), Roller Clamp (7), Male Luer Lock Adaptor (8), Male Luer Lock Anti-Free Flow Valve (9), Luer Lock Tip protector (10),



**P/N: A2-83072-N048**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Slide clamp (5)&(8), Two Needleless Injection Sites, one at the proximal side (6), one at the distal side (9), Roller Clamp (7), Male Luer Lock Adaptor (10), Male Luer Lock Anti-Free Flow Valve (11), Luer Lock Tip protector (12)

## Filter Sets




### P/N: A2-80070-F

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 0.22 micron IVEX filter (5), Male Luer Lock Adaptor (6), Luer Lock Tip protector (7)




### P/N: A2-80071-PF

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Needleless Injection Site (5) at proximal side, Roller clamp (6), 0.22 micron IVEX filter (7), Male Luer Lock Adaptor (8), Luer Lock Tip protector (9)



### P/N: A2-80071-DF


105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 0.22 micron IVEX filter (5), Slide clamp (6), Needleless Injection Site (7) at distal side, Male Luer Lock Adaptor (8), Luer Lock Tip protector (9)



### P/N: A2-80072-F

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (9), Roller Clamp (6), 0.22 micron IVEX filter (7), Slide clamp (8), Male Luer Lock Adaptor (10), Luer Lock Tip protector (11)

# Filter Sets




P/N: A2-80072-N248

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Slide clamp (5)&(9), Two Needless Injection Sites, one at the proximal side (6), one at the distal side (10), Roller Clamp (7), 0.22 micron IVEX filter (8), Male Luer Adaptor (11), Luer Lock Tip protector (12)

# Filter Sets



**P/N: A2-84071-D258**


105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 0.22 micron IVEX filter (5), Slide clamp (6), Needleless Injection Site (7) at distal side, Male Luer Lock Adaptor (8), Vented Cap with Filter (9)



**P/N: A2-84072-N248**


105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Slide clamp (5) & (9), Two Needleless Injection Sites, one at the proximal side (6), one at the distal side (10), Roller Clamp (7), 0.22 micron IVEX filter (8), Male Luer Lock Adaptor (11), Vented Cap with Filter (12)

## Secondary Set




**P/N: A2-80075**

40" (101cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller clamp (4), Male Luer Lock Adaptor (5), Luer Lock Tip protector (6), Bag Hanger (7)



**P/N: A2-84075-N056**

40" (101cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller clamp (4), Male Luer Lock Adaptor (5), Vented Cap with Filter (6), Bag Hanger (7)



**P/N: A2-80071-DFAC**

80" (203cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free  
 Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 0.22 micron IVEX filter (5), Slide clamp (6), Needless Injection Site (7) at distal side, Male Luer Lock Adaptor (8), Luer Cap with Filter (9)


## Filter Sets with Anti-Free-Flow Valve



**P/N: A2-83070-F**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller clamp (4), 0.22 micron IVEX filter (5), Male Luer Lock Adaptor (6), Male Luer Lock Anti-Free-Flow Valve (7), Luer Lock Tip protector (8)



**P/N: A2-83071-PF**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Needleless Injection Site (5) at proximal side, Roller clamp (6), 0.22 micron IVEX filter (7), Male Luer Lock Adaptor (8), Male Luer Lock Anti-Free-Flow Valve (9), Luer Lock Tip protector (10)



**P/N: A2-83071-DF**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 0.22 micron IVEX filter (5), Slide clamp (6), Needleless Injection Site (7), Male Luer Lock Adaptor (8), at distal side, Male Luer Lock Anti-Free-Flow Valve (9), Luer Lock Tip protector (10)




**P/N: A2-83072-F**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (9), Roller Clamp (6), 0.22 micron IVEX filter (7), Slide clamp (8), Male Luer Lock Adaptor (10), Male Luer Lock Anti-Free-Flow Valve (11), Luer Lock Tip protector (12)

# Filter Sets with Anti-Free-Flow Valve




**P/N: A2-83072-N248**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Slide clamp (5)&(9), Two Needleless Injection Sites, one at the proximal side (6), one at the distal side (10), Roller Clamp (7), 0.22 micron IVEX filter (8), Male Luer Lock Adaptor (11), Male Luer Lock Anti-Free-Flow Valve (12), Luer Lock Tip protector (13)


## Filter Set



**P/N: A2-80070-F-120**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 1.2 micron Adult I.V. filter (5), Male Luer Lock Adaptor (6), Luer Lock Tip protector (7)



**P/N: A2-80071-PF-120**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free


Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Needleless Injection Site (5) at proximal side, Roller clamp (6), 1.2 micron Adult I.V. filter (7), Male Luer Lock Adaptor (8), Luer Lock Tip protector (9)



**P/N: A2-80071-DF-120**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Roller Clamp (4), 1.2 Micron Adult I.V. filter (5), Slide clamp (6), Needleless Injection Site (7) at distal side, Male Luer Lock Adaptor (8), Luer Lock Tip protector (9)




**P/N: A2-80072-F-120**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (9), Roller Clamp (6), 1.2 micron Adult I.V. filter (7), Slide clamp (8) Male Luer Lock Adaptor (10), Luer Lock Tip protector (11)

# Filter Sets



**P/N: A2-83072-F-120**

105" (266cm) Long, Fluid path is sterile, non-Pyrogenic, non-DEHP, Latex-free

Spike tip protector (1), Universal spike (2), Drip chamber (3), Back-Check Valve (4), Two Needleless Injection Sites, one at the proximal side (5), one at the distal side (9), Roller Clamp (6), 1.2 micron Adult I.V. Filter (7), Slide clamp (8), Male Luer Lock Adaptor (10), Male Luer Lock Anti-Free-Flow Valve (11), Luer Lock Tip protector (12)